

Sonoma County Public Health Laboratory

3313 Chanate Road
Santa Rosa, CA 95404
(707) 565-4711

Quantiferon Gold Plus Testing

Patient Information

Last Name _____ First Name _____

Age _____ Date of Birth _____ Gender: M F

Medical Record Number _____ Submitter Specimen Number _____

Check one:

- Tuberculosis exposure Occupational Screening

Specimen type and/or specimen source

≥ 5 ml lithium heparin tube. Stable at room temperature for only 3 hours. Promptly refrigerate at 4-8°C and send on cold packs to the laboratory within 36 hours of the date collected.

Date Collected _____ **Time Collected** _____

Physician's Name _____

Submitting Facility's name and complete mailing address

